

Carnation Ambulatory Monitor (CAM) Owner Instructions:

- The Carnation Ambulatory Monitor (CAM) is a rugged, very lightweight, comfortable device and in very few cases will the test be ruined. If we have a reason to repeat it, we can do so right away without any harm to the patient.
- The CAM is *water resistant* (not water proof), so rain/splash of water is not a concern.
- Do not let the patient swim or receive a bath during the recording period.
- Use the included diary to document any events (collapse, weakness, 'dazed' appearance, etc.) – always noting the time the event occurred.
 - **The included diary MUST be returned in the shipping box.**
- If the bandage slips, please try to adjust them to keep the device covered and protected. If this is impossible, please return to your trusted veterinary office for an assessment and to determine if early removal is required. If the study is incomplete, we can repeat the test if necessary. There may be an additional charge for repeat tests.
 - **We strongly recommend a Thunder Shirt (or similar) to better protect the device from any scratching or other damage.** These are available at all local pet stores and are excellent at protecting the device, ***preventing repeat testing.***
- In the event of detachment/problems with the monitor:
 - Please contact us via info@olympicvetcardio.com to discuss the situation.
 - We may still have enough data to establish a definitive diagnosis.
 - If the Recorder (small part with a silver button on it) detaches from the Battrode (part adhered to the skin), please return to your trusted veterinarian as soon as possible.
 - Please note the *time the Recorder was detached* (or found detached) and the time it was *reattached* in the diary. This does not void the study.
 - To reattach the Recorder, simply clip it back in see the LED light flash - potentially secure with medical tape.

PLEASE NOTE

**Due to the detailed analysis performed on each recording,
it may take up to 10-14 days to obtain results.
We will send a report as soon as it is available.**

Length of Recording: 48 hour / 7 Day

REMOVAL DATE: _____ / TIME: _____

Pet Owners – Please follow your veterinarian's instructions and return to your veterinarian's clinic at the above date and time for removal OR follow the included instructions for removal.